

ELECTRO CAUTERY



Electro Cauterization is the process of destroying tissue using heat conduction from a metal probe heated by electric current. The procedure is used to stop bleeding from small vessels (larger vessels being ligated) or for cutting through soft tissue and Cauterization was used to stop heavy bleeding.

Special Features:

- Offers good choice of monopolar and bipolar current.
- Portable model with briefcase for safe handling.
- Intensity variable accurately (50 to 100 watts output)
- Ideal for minor and instant coagulation, dissection and fulguration.
- Useful for Dermatologists, ENT Surgeons, General surgeons and Cosmetologists.
- This method is preferred by surgeons by removal of warts Acne, Keloid, Corns, cysts, chonicrood, small tumours, conductivities etc.
- Excellent cosmetic results obtained with little or no side effects.
- Abnormalities are removed bloodlessly.

Standard Accessories :



Note: Unit is available with optional RF mode.

M/s. Healkeren Medtronics Private Limited

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

www.healkeren.com