

Infant Radiant Warmer with detachable baby trolley

Model : Heal - RWDT

III CE  
ISO 9001:2008  
ISO 13485:2003



Healiker<sup>TM</sup>

## Infant Radiant Warmer with detachable baby trolley : Heal - RWDT

Infant Radiant Warmer with detachable baby trolley (Heal - RWDT) is designed innovatively to meet neonatal care in well run and busy NICU units. The Detachable baby trolley can be used along with the warmer and it can be used separately as a baby cradle. Microprocessor based Temperature controller system provided to control the baby body temperature and it has all Fault Indications with visual alarms.

### TECHNICAL SPECIFICATIONS

#### SILENT FEATURES :

- Self Test System
- Auto/Manual selection modes
- Air/Skin servo controlled system
- Double Digital Display for ACTUAL & SET temperature
- °C/°F settings
- Heater power indicator
- Feather touch keys for easy operation
- Key pad lock facility
- AIR/SKIN/MANUAL Modes
- Power failure alarm -Fault indication alarms for 1°C low/high, probe failures, high temperature cutoff, system failure, heater failure etc.,

#### PHYSICAL PARAMETERS :

Baby Tray	
Length	: 700mm
Width	: 500mm
Mattress Thickness	: 25mm
Over all	
Length	: 900mm
Width	: 800mm
Height	: 1900mm
Working level from ground	: 900mm
Weight Approx.	: 40Kg
Castors	: 75mm with 2 brake wheels and 2 non-brake wheels

#### BABY & AIR PROBE :

Baby & Air probe is made of highly sensitive sensor and is interchangeable, field calibration taken care of by software. The probe is made up of silicon sheathed cable, which allows easy cleaning.

#### RECOMMENDED ENVIRONMENTAL SPECIFICATIONS :

Operating temperature range	: 20°C to 33°C
Storage temperature range	: -25°C to 60°C
Operating humidity range	: 0 - 100% RH

#### OBSERVATION LAMP :

LED based observation light focus on any area on the bed.

#### HEATER :

Single element CE approved ceramic heater placed in a parabolic reflector.

#### SKIN / AIR TEMPERATURE PROBE :

Range	: 28°C to 38°C
Accuracy	: +/-0.2°C within range
Resolution	: 0.1°C
Probe interchangeability.	: +/-0.2°C (Servo)
Heater output controllable from 0 to 100% in Manual mode (continuously variable)	

#### BABY TROLLEY :

The Detachable baby trolley which consists of large baby tray & is made of mild steel frame with three side acrylic collapsible walls. Trendelenburg & Fowler system provided.

#### ELECTRICAL SPECIFICATIONS :

Electrical Supply	: 230 V $\pm$ 10%, 50 Hz
Power consumption with Heater at max. Output	: 650W

#### OTHER FACILITIES :

- I V stand provided
- Monitor tray provided

## M/s. Healkeren Medtronics Private Limited

No.4, 2nd Cross, Lakshmipura Main Road, Krishnappa Industrial Estate, Abbigere, Bengaluru - 560 090, Karnataka, India.

Mobile : 93101 18774 / 80882 84888 / 80882 84898, Email : sales@healkeren.com / healkeren@gmail.com / info@healkeren.com

[www.healkeren.com](http://www.healkeren.com)